

# UltraTRAK™

BLOOD GLUCOSE MONITORING SYSTEM

## Quick Reference Card



STEP 1. INSERT A TEST STRIP



STEP 2. MATCH THE CODE NUMBER



STEP 3. APPLY A DROP OF BLOOD



STEP 4. RESULTS IN 7 SECONDS


Accurate Results You Can Trust

# UltraTRAK™

BLOOD GLUCOSE MONITORING SYSTEM

## Quick Reference Card

### Problem Solving Guide

Message	Meaning	Action
E-b	The  symbol appears with this message on the display by itself. This means that the battery will not provide enough power for a test.	Replace the battery at once.
E-U	Error message could be caused by a used test strip, or a temporary or permanent electronics problem occurred.	Repeat the test with a new test strip. If the error message appears again, please call your local agent for help.
E-t	The temperature of the environment, meter, or test strip was below or above the system operation range. You can not perform a test until the meter and test strip reach a temperature within the operation range of 50~104°F (10~40°C).	Repeat the test after the meter and test strip have reached a temperature within the operation range.
E-0	Error message indicates that there may be a problem with the meter or code strip.	Repeat the test with a new code strip. If the problem persists, please call the Customer Care Line for help.
E-A	Error message indicates that there is a problem with the meter.	Review the instructions and try again with a new test strip. If the problem persists, please call your local agent for help.
E-E		
E-9	Error message indicates that you may remove the strip after applying blood to the absorbent hole.	Review the instructions and try again with a new test strip.

Accurate Results You Can Trust