

# Live your Life



Ian Brenner, Age 54  
Underwater Photographer and Diabetes Patient

## Managing Your Diabetes Just Got Easier, So You Can Live Your Life to the Fullest

### UltraTRAK™ PRO

Blood Glucose Monitoring System


*Rapid · Accurate · Reliable*

No Coding Necessary  
For the Simplest Testing & Accuracy

Strip Ejection  
For the Safest Handling

Uses AAA Batteries  
Easy to Replace

Alternate Site Testing  
More Options, Less Pain



Rapid Results  
In Just 4 to 7 Seconds

Tiny 0.7ul Blood Sample  
Virtually Painless Experience

Four Reminder Alarms  
Alerts You Can Count On

5000 Test Determinants  
Extensive Battery Life

UltraTRAK™ PRO is the most technologically advanced systems on the market. And with competitive pricing, unparalleled customer service, and a five year warranty, we make it easier for you to experience effective diabetes management.



For more information, please visit [www.vertexdiagnostics.com](http://www.vertexdiagnostics.com) or call 888.90.VERTEX

# Ultra**TRAK**<sup>TM</sup> PRO

Ultra**TRAK**<sup>TM</sup> PRO delivers the most advanced features for effective diabetes management for both point-of-care and self-testing



- Strip Ejection Capability
- No Coding Necessary
- Rapid Results in 4 to 7 Seconds
- Tiny 0.7ul Blood Sample Size
- 5000 Test Determinants
- 450 Test Memory with Date/Time
- Uses AAA Batteries
- Four Reminder Alarms
- 7.14.21.28.60.90 Day Averaging
- Stores Control Solution Testing Results
- Ketone Warning and Alternate Site Testing
- Sleek Design, Simple One Button Operation
- Five Year Warranty

Product Name	Product Number
Ultra <b>TRAK</b> PRO Meter	UTPBGM
Ultra <b>TRAK</b> PRO Test Strips 50/box	UTPBGS
Ultra <b>TRAK</b> PRO Control Solution (Low and High Levels)	UTPCS




Manufactured by Vertex Diagnostics. Please visit [www.vertexdiagnostics.com](http://www.vertexdiagnostics.com) for more information.

**Easy and Effective Diabetic Testing So You Can  
*Live Your Life***


# Ultra *TRAK*<sup>TM</sup> PRO

## Blood Sample Size Comparison



Blood drops not actual size

## Results Timing Comparison



## Alternate Site Testing

Ultra *TRAK*<sup>TM</sup> PRO allows for testing on the palms, forearms, upper arms, calves, and thighs with equivalent results to fingertip testing.

